



TEMPORARY ABUTMENT/ OPEN TRAY COPING

Platform	\bigcirc		Coping screw
Ø4.1	IPD/AA-PR-06	IPD/AA-PR-07	IPD/AA-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using





OPEN AND CLOSED TRAY COPING

Platform		Reference
ø3.5 Tiny®	Open tray coping	IPD/GA-CN-00
4/1	Open tray coping	IPD/AA-CR-00
ø4.1	Closed tray coping	IPD/AA-CR-01

Classic design.



HEALING ABUTMENT

Platform	Height (mm)	Reference
ø3.5 Tiny®	4.0	IPD/GA-DN-04
	2.0	IPD/AA-DR-02
	3.0	IPD/AA-DR-03
φ4.1 ———	4.0	IPD/AA-DR-04
	5.0	IPD/AA-DR-05



STRAIGHT AND ANGLED ABUTMENT Straight, 15° & 25°

Platform	Height (mm)	Angulation	\Box	$\square\square\square$
ø3.5 Tiny®		Straight	IPD/GA-PN-00	-
		Straight	IPD/AA-PR-00	IPD/AA-PR-01
	1.0	15°	IPD/AA-PR-02	-
ø 4.1	1.0	25°	IPD/AA-PR-03	-
	2.0	15°	IPD/AA-PR-04	-
	2.0	25°	IPD/AA-PR-05	-
45.5		Straight	IPD/GA-PW-00	-
Ø5.5		Straight	IPD/GA-PW-01	-

The IPD/GA-PW-01 model has a broader shoulder than the IPD/GA-PW-00. Using the low profile IPD/AA-TR-01 screw for IPD/AA-PR-03 is recommended.





CASTABLE

Platform	\bigcirc	
ø4.1	IPD/AA-HR-00	IPD/AA-RR-00
ø5.5	IPD/GA-HW-00	IPD/GA-RW-00

Compatibility

BTI®

GB Series **INTERNA UNIVERSAL®**

[Universal® ø4.1 / Wide ø5.5]







ANALOG / DIGITAL ANALOG

COMPATIBILITY

Platform	Reference	Digital reference
ø4.1	IPD/GB-AR-00	IPD/GB-AR-00/3D
ø5.5	IPD/GB-AW-00	-
Lower Screw	-	IPD/KA-TA-00
Transverse Screw	-	IPD/KA-TA-01





SCREW

Platform	Screwdriver	Seat	Reference	TiN Reference
ø4.1 ø5.5	1.20 Hex	Flat	IPD/GB-TR-00	IPD/GB-TR-00/TIN
HE	Ex 1.20 Lengths 18	/25/32 mm. Please se	ee page 224 for more information.	





TPA SCREW FOR ANGLED PROSTHESIS

Platform	Screwdriver	Seat	Reference
ø4.1 ø5.5	TPA	Flat	IPD/GB-TR-50
	Langths 18	25/22 mm Plaasa saa naga 226	for more information



INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Platform	Height (mm)	\Box				
4/4	10	IPD/GB-SR-00	IPD/GB-SR-01			
Ø4.1	15	IPD/GB-SR-02	-			





Captive screw included. 1.20 Hex Screwdriver.



Platform	Height (mm)	\bigcirc	
	0.5	IPD/GB-IR-00/3D	IPD/GB-IR-01/3D
ø4.1	1.5	IPD/GB-IR-02/3D	IPD/GB-IR-03/3D
	2.5	IPD/GB-IR-04/3D	IPD/GB-IR-05/3D



High version of the non-engaging Custom Interface subject to availability. Please see page 188 for the (pd) Custom Interface System® accessories.

STRAIGHT AND ANGLED MULTI-UNIT ABUTMENT

Straight, 17° & 30°

Height (mm)	Straight	17º	30°
1.0	IPD/GB-MR-01	-	-
2.0	IPD/GB-MR-02	-	-
3.0	IPD/GB-MR-03	IPD/GB-MR-13	-
4.0	IPD/GB-MR-04	IPD/GB-MR-14	IPD/GB-MR-34
5.0	IPD/GB-MR-05	-	IPD/GB-MR-35
	1.0 2.0 3.0 4.0	1.0 IPD/GB-MR-01 2.0 IPD/GB-MR-02 3.0 IPD/GB-MR-03 4.0 IPD/GB-MR-04	1.0 IPD/GB-MR-01 - 2.0 IPD/GB-MR-02 - 3.0 IPD/GB-MR-03 IPD/GB-MR-13 4.0 IPD/GB-MR-04 IPD/GB-MR-14





Our straight multi-units can be used for single and multiple restorations. Please see page 198 for the Multi-unit system accessories.

STRAIGHT AND ANGLED CASTABLE ABUTMENT WITH A CO-Cr BASE

Straight, 150 & 250

Platform	Angulation	\bigcirc	
	Straight	IPD/GB-BR-00	IPD/GB-BR-01
Ø4.1	15 ⁰	IPD/GB-BR-10	IPD/GB-BR-11
	25°	IPD/GB-BR-20	IPD/GB-BR-21





All angled castable abutments include a TPA screw.

Co-Cr 3D BASE / SCAN BODY

Platform		\Box		Scan Body	
ø4.1	Straight / Angled	IPD/GB-BR-00/3D	IPD/GB-BR-01/3D	IPD/SC-AN-B1	



PSD OVERDENTURE ABUTMENT

Platform	Height (mm)	Reference
	1.0	IPD/GB-LR-01
Ø4.1	3.0	IPD/GB-LR-03
	5.0	IPD/GB-LR-05



Please see page 208 for the PSD system accessories.

GB Series INTERNA UNIVERSAL®
[Universal® \$\psi4.1 \frac{1}{2}\$ Wide \$\psi5.5]



TEMPORARY ABUTMENT/ OPEN TRAY COPING

Platform	\bigcirc		Coping screw
ø4.1	IPD/GB-PR-06	IPD/GB-PR-07	IPD/GB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

Platform	Height (mm)	Reference
	1.0	IPD/GB-DR-01
φ4.1	3.0	IPD/GB-DR-03
	5.0	IPD/GB-DR-05



CASTABLE

Platform	\Box	
ø4.1	IPD/GB-HR-00	IPD/GB-RR-00
ø5.5	-	IPD/GB-RW-00



Compatibility

BTI®

GC Series MULTI-IM®

[ø4.1]



