





TEMPORARY ABUTMENT/ OPEN TRAY COPING

Platform	\bigcirc		Coping screw
ø3.4	IPD/BA-PN-06	IPD/BA-PN-07	IPD/AA-TR-00/C
ø4.1	IPD/AA-PR-06	IPD/AA-PR-07	IPD/AA-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



OPEN AND CLOSED TRAY COPING

	Reference
Open tray coping	IPD/AA-CR-00
Closed tray coping	IPD/AA-CR-01
Open tray coping	IPD/BA-CW-00
	Closed tray coping

Classic design.



HEALING ABUTMENT

Platform	Height (mm)	Reference
	2.0	IPD/BA-DN-02
ø3.4	3.0	IPD/BA-DN-03
	4.0	IPD/BA-DN-04
	2.0	IPD/AA-DR-02
4/1	3.0	IPD/AA-DR-03
ø4.1	4.0	IPD/AA-DR-04
	5.0	IPD/AA-DR-05
	2.0	IPD/BA-DW-02
ø5.0	3.0	IPD/BA-DW-03
	4.0	IPD/BA-DW-04



STRAIGHT AND ANGLED ABUTMENT

Straight, 15° & 25°

Platform	Height (mm)	Angulation	\Box	$\square\square$
		Straight	IPD/AA-PR-00	IPD/AA-PR-01
	1.0	15º	IPD/AA-PR-02	-
Ø4.1	1.0	25 ⁰	IPD/AA-PR-03	-
	2.0	15°	IPD/AA-PR-04	-
	2.0	25 ⁰	IPD/AA-PR-05 -	-

Using the low profile IPD/AA-TR-01 screw for IPD/AA-PR-03 is recommended.



CASTABLE

Platform	\bigcirc	ada
ø4.1	IPD/AA-HR-00	IPD/AA-RR-00
ø5.0	IPD/BA-HW-00	IPD/BA-RW-00



Compatibility

BIOMET® 3i®

BB Series **CERTAIN®**

[ø3.4 / ø4.1 / ø5.0]







CATALOG 2021









ANALOG / DIGITAL ANALOG

erence	Digital reference
/BB-AN-00	IPD/BB-AN-00/3D
/BB-AR-00	IPD/BB-AR-00/3D
/BB-AW-00	IPD/BB-AW-00/3D
	IPD/KA-TA-00
	IPD/KA-TA-01





SCREW

Platform	Screwdriver	Seat	Reference	TiN Reference
ø3.4 ø4.1 ø5.0	1.20 Hex	Flat	IPD/BB-TR-01	IPD/BB-TR-01/TIN
ø4.1 ø5.0	1.20 Hex	Conical	IPD/BB-TR-02	-

IPD/BB-TR-02 is a specific screw for models IPD/BB-RR-00 and IPD/BB-RW-00

HEX 1.20 Lengths 18/25/32 mm. Please see page 224 for more information.



TPA SCREW FOR ANGLED PROSTHESIS

Platform	Screwdriver	Seat	Reference
ø3.4 ø4.1 ø5.0	TPA	Flat	IPD/BB-TR-51



Lengths 18/25/32 mm. Please see page 226 for more information.



INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Platform	Height (mm)	\bigcirc	
42 /	10	IPD/BB-SN-00	IPD/BB-SN-01
Ø3.4	15	IPD/BB-SN-02	-
4/1	10	IPD/BB-SR-00	IPD/BB-SR-01
ø4.1	15	IPD/BB-SR-02	-
dE O	10	IPD/BB-SW-00	IPD/BB-SW-01
ø5.0	15	IPD/BB-SW-02	-

Captive screw included. 1.20 Hex Screwdriver.





INTERFACE (ipd) Custom Interface

	Dy Stein		
Platform	Height (mm)	\bigcirc	
	0.5	IPD/BB-IN-00/3D	IPD/BB-IN-01/3D
Ø3.4	1.5	IPD/BB-IN-02/3D	IPD/BB-IN-03/3D
_	2.5	IPD/BB-IN-04/3D	IPD/BB-IN-05/3D
4/1	0.5	IPD/BB-IR-00/3D	IPD/BB-IR-01/3D
Ø4.1 –	1.5	IPD/BB-IR-02/3D	IPD/BB-IR-03/3D
45.0	0.5	IPD/BB-IW-00/3D	IPD/BB-IW-01/3D
φ5.0 –	1.5	IPD/BB-IW-02/3D	IPD/BB-IW-03/3D

High version of the non-engaging Custom Interface subject to availability. Please see page 188 for the (ipd) Custom Interface System® accessories.



STRAIGHT AND ANGLED MULTI-UNIT ABUTMENT

Platform	Height (mm)	Straight	17º	30°
	1.0	IPD/BB-MR-01	-	-
	2.0	IPD/BB-MR-02	IPD/BB-MR-12	-
ø4.1	3.0	IPD/BB-MR-03	IPD/BB-MR-13	-
	4.0	IPD/BB-MR-04	-	IPD/BB-MR-34
	5.0	IPD/BB-MR-05	=	-

Our straight multi-units can be used for single and multiple restorations. Please see page 198 for the Multi-unit system accessories.



STRAIGHT AND ANGLED CASTABLE ABUTMENT WITH A Co-Cr BASE

Straight, 15° &	25°		
Platform	Angulation	\bigcirc	
	Straight	IPD/BB-BN-00	IPD/BB-BN-01
Ø3.4	15°	IPD/BB-BN-10	IPD/BB-BN-11
	25°	IPD/BB-BN-20	IPD/BB-BN-21
	Straight	IPD/BB-BR-00	IPD/BB-BR-01
Ø4.1	15°	IPD/BB-BR-10	IPD/BB-BR-11
	25°	IPD/BB-BR-20	IPD/BB-BR-21

All angled castable abutments include a TPA screw.



37









Co-Cr 3D BASE / SCAN BODY

Ø3.4 Straight IPD/BB-BN-00/3D IPD/BB-BN-01/3D IPD/SC	ody
W3.4	-AN-B1
Angled IPD/BB-BN-A0/3D IPD/BB-BN-A1/3D IPD/SC	-AN-B1
φ4.1 Straight / Angled IPD/BB-BR-00/3D IPD/BB-BR-01/3D IPD/SC	-AN-B1



PSD OVERDENTURE ABUTMENT

Platform	Height (mm)	Reference
	1.0	IPD/BB-LN-01
ø3.4	3.0	IPD/BB-LN-03
	5.0	IPD/BB-LN-05
	1.0	IPD/BB-LR-01
ø4.1	3.0	IPD/BB-LR-03
	5.0	IPD/BB-LR-05



Please see page 208 for the PSD system accessories.

TEMPORARY ABUTMENT/ OPEN TRAY COPING

Platform	\Box		Coping screw
ø3.4	IPD/BB-PN-06	IPD/BB-PN-07	IPD/BB-TR-01/C
ø4.1	IPD/BB-PR-06	IPD/BB-PR-07	IPD/BB-TR-01/C



The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.

OPEN AND CLOSED TRAY COPING

Platform		Reference
ø5.0	Open tray coping	IPD/BB-CW-00



HEALING ABUTMENT

Platform	Height (mm)	Reference
	2.0	IPD/BB-DN-02
Ø3.4	3.0	IPD/BB-DN-03
	4.0	IPD/BB-DN-04
	2.0	IPD/BB-DR-02
ø4.1	3.0	IPD/BB-DR-03
	4.0	IPD/BB-DR-04
	2.0	IPD/BB-DW-02
ø5.0	3.0	IPD/BB-DW-03
	4.0	IPD/BB-DW-04



STRAIGHT ABUTMENT

Platform	Reference
Ø3.4	IPD/BB-PN-00
ø4.1	IPD/BB-PR-00
ø5.0	IPD/BB-PW-00





CASTABLE

Platform	\bigcirc	
ø4.1	IPD/BB-HR-00	IPD/BB-RR-00
ø5.0	IPD/BB-HW-00	IPD/BB-RW-00

Engaging castables should be used with the IPD/BB-TR-02 screw.

