




Compatibility



**CAMLOG®**

**JB Serie CONELOG®**

[ $\phi 3.3$  /  $\phi 3.8$  /  $\phi 4.3$ ]





## INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Platform	Height (mm)		
$\phi 3.3$	10	IPD/JB-ST-00	IPD/JB-ST-01
	15	IPD/JB-ST-02	-
$\phi 3.8$	10	IPD/JB-SN-00	IPD/JB-SN-01
	15	IPD/JB-SN-02	-
$\phi 4.3$	10	IPD/JB-SR-00	IPD/JB-SR-01
	15	IPD/JB-SR-02	-

Captive screw included. 1.20 Hex Screwdriver.



## Ti-BASE Custom Ti-Base System®

Platform	Height (mm)		
$\phi 3.3$	0.5	IPD/JB-IT-00/3D	IPD/JB-IT-01/3D
	1.5	IPD/JB-IT-02/3D	IPD/JB-IT-03/3D
$\phi 3.8$	0.5	IPD/JB-IN-00/3D	IPD/JB-IN-01/3D
	1.5	IPD/JB-IN-02/3D	IPD/JB-IN-03/3D
$\phi 4.3$	0.5	IPD/JB-IR-00/3D	IPD/JB-IR-01/3D
	1.5	IPD/JB-IR-02/3D	IPD/JB-IR-03/3D

Please see [page 147](#) for the  Custom Ti-Base System® accessories.



## SCREW

Platform	Screwdriver	Seat	Reference	TiN Reference
$\phi 3.3$ $\phi 3.8$ $\phi 4.3$	Hex 1.25	Cónico	IPD/JB-TR-00 IPD/JB-TR-01	IPD/JB-TR-00/TiN

The reference IPD/JB-TR-01 is an abutment extractor screw.

The IPD/JB-TR-01 Abutment Extractor Screw cannot be used if angled systems are used with the ConeLog® System.

 Lengths 18/25/32 mm. Please see [page 184](#) for more information.



## ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
$\phi 3.3$	IPD/JB-AT-00	IPD/JB-AT-00/3D
$\phi 3.8$	IPD/JB-AN-00	IPD/JB-AN-00/3D
$\phi 4.3$	IPD/JB-AR-00	IPD/JB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





STRAIGHT AND ANGLED MULTI-UNIT ABUTMENT  
Straight, 17° & 30°

Platform	Height (mm)	Straight	17°	30°
ø3.8	1.5	IPD/JB-MN-01	-	-
	2.5	IPD/JB-MN-02	IPD/JB-MN-12	-
	3.5	-	-	IPD/JB-MN-33
ø4.3	1.5	IPD/JB-MR-01	-	-
	2.5	IPD/JB-MR-02	IPD/JB-MR-12	-
	3.5	-	-	IPD/JB-MR-33

Our straight Multi-units can be used for single and multiple restorations.  
Please see [page 157](#) for the Multi-unit system accessories.



TEMPORARY ABUTMENT / OPEN TRAY COPING

Platform			Coping screw
ø3.3	IPD/JB-PT-06	IPD/JB-PT-07	IPD/JB-TR-00/C
ø3.8	IPD/JB-PN-06	IPD/JB-PN-07	
ø4.3	IPD/JB-PR-06	IPD/JB-PR-07	

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



PSD OVERDENTURE ABUTMENT

Platform	Height (mm)	Reference
ø3.8	2.0	IPD/JB-LN-02
	3.0	IPD/JB-LN-03
	4.0	IPD/JB-LN-04
ø4.3	2.0	IPD/JB-LR-02
	3.0	IPD/JB-LR-03
	4.0	IPD/JB-LR-04

Please see [page 169](#) for the PSD system accessories.



HEALING ABUTMENT

Platform	Height (mm)	Reference
ø3.3	2.0	IPD/JB-DT-02
	4.0	IPD/JB-DT-04
	6.0	IPD/JB-DT-06
ø3.8	2.0	IPD/JB-DN-02
	4.0	IPD/JB-DN-04
	6.0	IPD/JB-DN-06
ø4.3	2.0	IPD/JB-DR-02
	4.0	IPD/JB-DR-04
	6.0	IPD/JB-DR-06

