

# Compatibility STRAUMANN® DD Serie SRA®

[ø3,5 / ø4,6]



### INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

| Platform | Height (mm) | $\bigcirc$ |              |
|----------|-------------|------------|--------------|
| ø3.5     | 8.5         | -          | IPD/DD-SN-01 |
| ø4.6     | 8.5         | -          | IPD/DD-SR-01 |

Captive screw included. 1.20 Hex Screwdriver.





# Ti-BASE (ipd) Custom Interface System\*

| Platform | Height (mm) | $\Box$ |                 |
|----------|-------------|--------|-----------------|
| ø3.5     | 0.2         | -      | IPD/DD-IN-01/3D |
| Ø4,6     | 0.2         | -      | IPD/DD-IR-01/3D |



#### **SCREW**

| Platform    | Screwdriver | Seat    | Reference    | TiN reference    |
|-------------|-------------|---------|--------------|------------------|
| ø3.5 / ø4.6 | -           | Conical | IPD/DD-TR-00 | IPD/DD-TR-00/TIN |
|             | Torx        |         | -            | IPD/DD-TR-50     |





# **ANALOG / DIGITAL ANALOG**

| Platform     | Reference    | Digital Reference |
|--------------|--------------|-------------------|
| ø3.5         | IPD/DD-AN-00 | IPD/DD-AN-00/3D   |
| ø4.6         | IPD/DD-AR-00 | IPD/DD-AR-00/3D   |
| Bottom screw | -            | IPD/KA-TA-00      |
| Side screw   | -            | IPD/KA-TA-01      |





# **TEMPORARY ABUTMENT**

| Platform | $\Box$       |              | Copying screw     |
|----------|--------------|--------------|-------------------|
| ø3.5     | IPD/DD-PN-06 | IPD/DD-PN-07 | - IPD/DD-TR-00/C  |
| Ø4.6     | IPD/DD-PR-06 | IPD/DD-PR-07 | = 1FD/DD-11C-00/C |

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.





# **HEALING ABUTMENT**

| Platform | Height (mm) | Reference    |
|----------|-------------|--------------|
| ø3.5     | 4.5         | IPD/DD-DN-04 |
| ø4.6     | 4.5         | IPD/DD-DR-04 |



