


Compatibility  
**NEODENT®**  
**RB Series**  
**GRAND MORSE™ GM**  
**[GM]**






COMPATIBILITY  
**RB Series**  
[GM]

**GRAND MORSE™ GM**



COMPATIBILITY  
**RB Series**

**GRAND MORSE™ GM**  
[GM]




## ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
GM	IPD/RB-AR-00	IPD/RB-AR-00/3D
Lower Screw	-	IPD/KA-TA-00
Transverse Screw	-	IPD/KA-TA-01



## SCREW


Platform	Screwdriver	Seat	Reference	TiN Reference
GM	Neo GM	Conical	IPD/RB-TR-00	IPD/RB-TR-00/TiN

 Lengths 18/25/32 mm. Please see page 225 for more information.





## TPA SCREW FOR ANGLED PROSTHESIS

Platform	Screwdriver	Seat	Reference
GM	TPA	Conical	IPD/RB-TR-50

 Lengths 18/25/32 mm. Please see page 226 for more information.





## INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Platform	Height (mm)		
GM	10	IPD/RB-SR-00	IPD/RB-SR-01
	15	IPD/RB-SR-02	-

Captive screw included. 1.20 Hex Screwdriver.



## INTERFACE

Platform	Height (mm)		
GM	1.0	IPD/RB-IR-00/3D	IPD/RB-IR-01/3D
	2.0	IPD/RB-IR-02/3D	IPD/RB-IR-03/3D
	3.0	IPD/RB-IR-04/3D	IPD/RB-IR-05/3D

Please see page 188 for the  Custom Interface System® accessories.



## STRAIGHT AND ANGLED MULTI-UNIT ABUTMENT

Straight, 17° & 30°

Platform	Height (mm)	Straight	17°	30°
GM	1.5	IPD/RB-MR-01	-	-
	2.5	IPD/RB-MR-02	IPD/RB-MR-12	-
	3.5	IPD/RB-MR-03	-	IPD/RB-MR-33

Our straight multi-units can be used for single and multiple restorations.  
Please see page 198 for the Multi-unit system accessories.




## STRAIGHT CASTABLE ABUTMENT WITH A Co-Cr BASE

Platform		
GM	IPD/RB-BR-00	IPD/RB-BR-01



## Co-Cr 3D BASE / SCAN BODY

Platform			Scan Body
GM	Straight / Angled	IPD/RB-BR-00/3D	IPD/RB-BR-01/3D



## PSD OVERDENTURE ABUTMENT

Platform	Height (mm)	Reference
GM	2.0	IPD/RB-LR-02
	3.0	IPD/RB-LR-03
	4.0	IPD/RB-LR-04

Please see page 208 for the PSD system accessories.



## TEMPORARY ABUTMENT/ OPEN TRAY COPING

Platform			Coping screw
GM	IPD/RB-PR-06	IPD/RB-PR-07	IPD/RB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



## HEALING ABUTMENT

Platform	Height (mm)	Reference
GM	1.5	IPD/RB-DR-01
	3.0	IPD/RB-DR-03
	5.0	IPD/RB-DR-05



 The registered brands in IPD 2004 are the property of IMPLANT PROTESIS DENTAL 2004 SL. The other brands mentioned are the property of their respective owners, and IMPLANT PROTESIS DENTAL 2004 SL has no relationship with either they or the aforementioned brands.