



Compatibility



STRAUMANN®

DA Serie TISSUE LEVEL

[ø4.8 RN / ø6.5 WN]





INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Platform	Height (mm)		
ø4.8 RN	10	IPD/DA-SR-00	IPD/DA-SR-01
	15	IPD/DA-SR-02	-
ø6.5 WN	10	IPD/DA-SW-00	IPD/DA-SW-01

Captive screw included. 1.20 Hex Screwdriver.

Ti-BASE Custom Ti-Base System®

Platform	Height (mm)		
ø4.8 RN	0.4	IPD/DA-IR-00/3D	IPD/DA-IR-01/3D
ø6.5 WN	0.4	IPD/DA-IW-02/3D	IPD/DA-IW-03/3D


Please see [page 147](#) for the  Custom Ti-Base System® accessories.





SCREW / TPA SCREW FOR ANGLED PROSTHESIS


Platform	Screwdriver	Seat	Reference	TiN Reference
ø4.8 RN ø6.5 WN Implant level	Torx	Conical	IPD/DA-TR-01	IPD/DA-TR-01/TiN
	Torx	Flat	IPD/DA-TR-02	-
	Torx	Conical	IPD/DA-TR-03	IPD/DA-TR-03/TiN
ø4.8 RN synOcta®	Torx	Conical	IPD/DA-TR-00	IPD/DA-TR-00/TiN
ø4.8 RN ø6.5 WN	TPA	Conical synOcta®	-	IPD/DA-TR-50
		Implant level conical	-	IPD/DA-TR-51

The IPD/DA-TR-03 screw has the same measurements as the original direct implant screw.

IPD/DA-TR-01 screw is a narrower headed screw to be used on  Custom Ti-Base System®

 Lengths 18/25/32 mm. Please see [page 185](#) for more information.

 Lengths 18/25/32 mm. Please see [page 186](#) for more information.



ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø4.8 RN	IPD/DA-AR-00	IPD/DA-AR-00/3D
	synOcta®	-
ø6.5 WN	IPD/DA-AW-00	IPD/DA-AW-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





MULTI-UNIT STRAIGHT ABUTMENT

Platform	Height (mm)	Straight
ø4.8 RN	0.9	IPD/DA-MR-01

Our straight Multi-units can be used for single and multiple restorations. Please see [page 157](#) for the Multi-unit system accessories.



TEMPORARY ABUTMENT / OPEN TRAY COPING

Platform				Coping screw
ø4.8 RN	synOcta®	IPD/DA-PR-06	IPD/DA-PR-07	IPD/DA-TR-01/C
		IPD/DA-PR-09	IPD/DA-PR-08	IPD/DA-TR-00/C
ø6.5 WN		IPD/DA-PW-06	IPD/DA-PW-07	IPD/DA-TR-01/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



PSD OVERDENTURE ABUTMENT

Platform	Height (mm)	Reference
ø4.8 RN	1.0	IPD/DA-LR-01
	2.0	IPD/DA-LR-02
	3.0	IPD/DA-LR-03

Please see [page 169](#) for the PSD system accessories.



HEALING ABUTMENT

Platform	Height (mm)	Reference
ø4.8 RN	0.0 Sealing cap	IPD/DA-DR-00
	1.5	IPD/DA-DR-01
	3.0	IPD/DA-DR-03
	4.5	IPD/DA-DR-04




STRAIGHT ABUTMENT SynOcta® and Straight

Platform	Reference
ø4.8 RN	IPD/DA-PR-00
	IPD/DA-PR-01
ø6.5 WN	IPD/DA-PW-00





CASTABLE

Platform			
ø4.8 RN	synOcta®	IPD/DA-HR-01	IPD/DA-RR-01
		IPD/DA-HR-00	IPD/DA-RR-00



STRAIGHT AND ANGLED CASTABLE ABUTMENT WITH A Co-Cr BASE Straight, 15° & 25°

Platform	Angulation		
ø4.8 RN	Straight	IPD/DA-BR-00	IPD/DA-BR-01
	15°	IPD/DA-BR-10	IPD/DA-BR-11
	25°	IPD/DA-BR-20	IPD/DA-BR-21
ø4.8 RN synOcta®	Straight	IPD/DA-BR-03	IPD/DA-BR-02
	15°	IPD/DA-BR-13	IPD/DA-BR-12
	25°	IPD/DA-BR-23	IPD/DA-BR-22
ø6.5 WN	Straight	IPD/DA-BW-00	IPD/DA-BW-01
	15°	IPD/DA-BW-10	IPD/DA-BW-11
	25°	IPD/DA-BW-20	IPD/DA-BW-21

Co-Cr Base for digital workflow	Scan Body
ø4.8 RN	IPD/SC-AN-B2
ø4.8 RN synOcta®	IPD/SC-AN-B2
ø6.5 WN	IPD/SC-AN-B2

All angled castable abutments include a TPA screw.

